

**Anti-DEC-205/CD205 antibody [NLDC145] {Alexa-fluor-546} (STJ16102647)**

STJ16102647

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-DEC-205/CD205 is suitable for use in FACS, Immunohistochemistry and Immunofluorescence research applications.
<b>Applications</b>	FACS/IHC-F/IF
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	NLDC145
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Alexa-fluor-546
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	
<b>Formulation</b>	TRIS, PBS, pH7 and antibiotics
<b>Isotype</b>	IgG2a
<b>Storage</b>	Store at -20°C for up to 1 year.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen</b>
<b>Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>