

Anti-MMP-12/HME antibody [603E6.22] {HRPO} (STJ16102648)

STJ16102648

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MMP-12/HME is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC-F/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	603E6.22
Concentration	0.5 mg/mL
Conjugation	HRPO
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence