

Anti-MMP-12/HME antibody [706F9.01] {HRPO} (STJ16102649)
STJ16102649**GENERAL INFORMATION**

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-MMP-12/HME is suitable for use in ELISA and Western Blot research applications. |
| Applications | ELISA/WB |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|----------------------------------|
| Clonality | Monoclonal |
| Clone ID | 706F9.01 |
| Concentration | 0.5 mg/mL |
| Conjugation | HRPO |
| Purification | Affinity Purified |
| Dilution Range | |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG1 |
| Storage Instruction | Store at -20°C for up to 1 year. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |