

Anti-MMP-12/HME antibody [701E4.03] {HRPO} (STJ16102651)
STJ16102651**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MMP-12/HME is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	701E4.03
Concentration	0.5 mg/mL
Conjugation	HRPO
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	