

Anti-MIP-3 Alpha/CCL20 antibody [206D9.05] {HRPO} (STJ16102656)
STJ16102656**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MIP-3 Alpha/CCL20 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	206D9.05
Concentration	0.5 mg/mL
Conjugation	HRPO
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence