

Anti-TLR8/CD288 antibody [307D3.01] {Alexa-fluor-647} (STJ16102658)

STJ16102658

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TLR8/CD288 is suitable for use in Western Blot and FACS research applications.
Applications	WB/FACS
Host/Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	307D3.01
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	170744
Gene Symbol	Tlr8
Uniprot ID	TLR8_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	