

Anti-AMACR antibody [13H4] (STJ180006)

STJ180006

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-AMACR is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	13H4
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:150
Formulation	Tris-HCl buffer containing stabilizing protein (BSA) and <0.1% ProClin
Isotype	IgG
Storage	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	23600
Gene Symbol	AMACR
Uniprot ID	AMACR_HUMAN
Immunogen	Human AMACR (P504S) polypeptide
Immunogen Region	
Specificity	Positive control: Prostate carcinoma
Immunogen Sequence	



Human prostate carcinoma stained anti-AMACR antibody using peroxidase-conjugate and DAB chromogen. Note cytoplasmic and luminal circumferential staining of carcinoma glands.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081