

Anti-BCA-225 antibody [Cu-18] (STJ180016)

STJ180016

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-BCA-225 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Cu-18
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:25-100
Formulation	Tris-HCl buffer containing stabilizing protein (BSA) and <0.1% ProClin
Isotype	IgG1/Kappa
Storage	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	A synthetic peptide
Immunogen Region	
Specificity	Positive control: Breast carcinoma, fallopian tube, or skin
Immunogen Sequence	



Human breast adenocarcinoma stained with anti-BCA-225 antibody using peroxidase-conjugate and DAB chromogen. Note the cytoplasmic staining of tumor cells.