

Anti-CD21 antibody (44-196aa) [ZM75] (STJ180316)

STJ180316

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD21 (44-196aa) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZM75
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:100-200
Formulation	Tris-HCl buffer containing stabilizing protein (BSA) and <0.1% ProClin
Isotype	IgG1/Kappa
Storage	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	1380
Gene Symbol	CR2
Uniprot ID	CR2_HUMAN
Immunogen	Recombinant human CR2 (CD21) protein fragment (around aa 44-196)
Immunogen Region	44-196aa
Specificity	Positive control: Tonsil, lymph node
Immunogen Sequence	



Human tonsil stained with anti-CD21 antibody using peroxidase-conjugate and DAB chromogen. Note the membrane staining of follicular dendritic cells in follicles.