

Anti-Erythropoietin antibody (28-162aa) [ZM135] (STJ180351)

STJ180351

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Erythropoietin (28-162aa) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZM135
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:100-200
Formulation	Tris-HCl buffer containing stabilizing protein (BSA) and <0.1% ProClin
Isotype	IgG1/Kappa
Storage	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	2056
Gene Symbol	EPO
Uniprot ID	EPO_HUMAN
Immunogen	Recombinant fragment of human EPO protein (around aa 28-162)
Immunogen	28-162aa
Region	
Specificity	Positive control: Kidney
Immunogen	
Sequence	



Human kidney stained with anti-Erythropoietin antibody using peroxidase-conjugate and DAB chromogen. Note cytoplasmic staining of tubular cells.