

Anti-Fascin antibody (Full length) [ZM192] (STJ180352)

STJ180352

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Fascin (Full length) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZM192
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:100-200
Formulation	Tris-HCl buffer containing stabilizing protein (BSA) and <0.1% ProClin
Isotype	IgG2a/Kappa
Storage	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	6624
Gene Symbol	FSCN1
Uniprot ID	FSCN1_HUMAN
Immunogen	Full length recombinant human FSCN1 protein
Immunogen Region	Full length
Specificity	Positive control: Hodgkin's lymphoma or spleen
Immunogen Sequence	



Human Hodgkin lymphoma stained with anti-Fascin antibody using peroxidase-conjugate and DAB chromogen. Note the cytoplasmic staining of Hodgkin NA's cells.