

Anti-ALK-p80 antibody (1058-1620aa) [ZR305] (STJ180426)

STJ180426

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-ALK-p80 (1058-1620aa) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZR305
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:100-200
Formulation	Tris-HCl buffer containing stabilizing protein (BSA) and <0.1% ProClin
Isotype	IgG
Storage	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human ALK-1 protein fragment around aa 1058-1620
Immunogen Region	1058-1620aa
Specificity	Positive control: Anaplastic large cell lymphoma or Lung adenocarcinoma with ALK translocation
Immunogen Sequence	



Human lung adenocarcinoma stained with anti-ALK-1 antibody using peroxidase-conjugate and DAB chromogen. Note cytoplasmic stain of adenocarcinoma cells.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081