

Anti-MUC-5AC antibody [ZR19] (STJ180524)

STJ180524

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-MUC-5AC is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZR19
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:100-200
Formulation	Tris-HCl buffer containing stabilizing protein (BSA) and <0.1% ProClin
Isotype	
Storage	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	4586
Gene Symbol	MUC5AC
Uniprot ID	MUC5A_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Formalin-fixed, paraffin-embedded human stomach stained with anti-MUC-5AC antibody using peroxidase-conjugate and DAB chromogen. Note the cytoplasmic staining of gastric body glands