

## Anti-Brachyury antibody (235-435aa) [ZR391] (STJ180550)

STJ180550

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit monoclonal antibody anti-Brachyury (235-435aa) is suitable for use in Immunohistochemistry research applications.
<b>Applications</b>	IHC-P
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	ZR391
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purified
<b>Dilution Range</b>	1:100-200
<b>Formulation</b>	Tris-HCl buffer containing stabilizing protein (BSA) and <0.1% ProClin
<b>Isotype</b>	IgG
<b>Storage</b>	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	6862
<b>Gene Symbol</b>	TBXT
<b>Uniprot ID</b>	TBXT_HUMAN
<b>Immunogen</b>	Recombinant fragment (around aa 235-435) of human TBXT protein
<b>Immunogen Region</b>	235-435aa
<b>Specificity</b>	Positive Control: Chordoma
<b>Immunogen Sequence</b>	



Formalin-fixed and paraffin-embedded human chordoma stained with anti-brachyury antibody using peroxidase-conjugate and DAB chromogen. Note the absence of nuclear staining of tumor cells.