

Anti-CD71 antibody [ZM136] (STJ180561)

STJ180561

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD71 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZM136
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:100-200
Formulation	Tris-HCl buffer containing stabilizing protein (BSA) and <0.1% ProClin
Isotype	IgG1/Kappa
Storage	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	7037
Gene Symbol	TFRC
Uniprot ID	TFR1_HUMAN
Immunogen	Recombinant fragment of human TFRC receptor
Region	
Specificity	Positive Control: Bone marrow
Immunogen Sequence	



Formalin-fixed, paraffin-embedded bone marrow stained with anti-CD71 antibody using peroxidase-conjugated and DAB chromogen. Note the cytoplasmic staining of erythroid precursors