

Anti-Cytokeratin 19 antibody (Full length) [ZR143] (STJ180563)

STJ180563

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Cytokeratin 19 (Full length) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZR143
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:100-200
Formulation	Tris-HCl buffer containing stabilizing protein (BSA) and <0.1% ProClin
Isotype	IgG
Storage	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	3880
Gene Symbol	KRT19
Uniprot ID	K1C19_HUMAN
Immunogen	Recombinant human full-length KRT19 protein
Immunogen Region	Full length
Specificity	Positive Control: Breast carcinoma
Immunogen Sequence	



Formalin-fixed, paraffin-embedded human pancreatic ductal carcinoma stained with anti-keratin 19 antibody using peroxidase-conjugate and DAB chromogen. Note the cytoplasmic staining of tumor cells