

Anti-PDCD1 antibody (1-200aa) [ZM357] (STJ180588)

STJ180588

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-PD-1 (1-200aa) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZM357
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:100-200
Formulation	Tris-HCl buffer containing stabilizing protein (BSA) and <0.1% ProClin
Isotype	IgG1/Kappa
Storage	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	29126
Gene Symbol	CD274
Uniprot ID	PD1L1_HUMAN
Immunogen	Recombinant fragment (around aa1-200) of human PDCD1 protein
Immunogen Region	1-200aa
Specificity	Positive Control: Lymph node
Immunogen Sequence	



Formalin-fixed, paraffin-embedded human tonsil stained with anti-PD-1 antibody using peroxidase-conjugate and DAB chromogen. Note cytoplasmic staining of T-cells