

## Anti-IgD antibody [ZR156] (STJ180611)

STJ180611

### GENERAL INFORMATION


|                          |                                                                                                        |
|--------------------------|--------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                     |
| <b>Short Description</b> | Rabbit monoclonal antibody anti-IgD is suitable for use in Immunohistochemistry research applications. |
| <b>Applications</b>      | IHC-P                                                                                                  |
| <b>Host/Source</b>       | Rabbit                                                                                                 |
| <b>Reactivity</b>        | Human                                                                                                  |

### PRODUCT PROPERTIES

|                       |                                                                                                                                                                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>      | Monoclonal                                                                                                                                                                                                                                         |
| <b>Clone ID</b>       | ZR156                                                                                                                                                                                                                                              |
| <b>Concentration</b>  |                                                                                                                                                                                                                                                    |
| <b>Conjugation</b>    | Unconjugated                                                                                                                                                                                                                                       |
| <b>Purification</b>   | Affinity purified                                                                                                                                                                                                                                  |
| <b>Dilution Range</b> | 1:100-200                                                                                                                                                                                                                                          |
| <b>Formulation</b>    | Tris-HCl buffer containing stabilizing protein (BSA) and <0.1% ProClin                                                                                                                                                                             |
| <b>Isotype</b>        |                                                                                                                                                                                                                                                    |
| <b>Storage</b>        | Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C. |
| <b>Instruction</b>    |                                                                                                                                                                                                                                                    |

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Immunogen Region**  
**Specificity**  
**Immunogen Sequence**



Formalin-fixed, paraffin-embedded human lymph node stained with anti-IgD antibody using peroxidase-conjugate and DAB chromogen. Note cytoplasmic staining of mantle B cells