

## Anti-BCL3 antibody (340-420 aa) (STJ190158)

STJ190158

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-B-cell lymphoma 3 protein (340-420 aa) is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
<b>Applications</b>	IHC/IF/WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Rat/Mouse/Goat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution Range</b>	WB 1:500-2000 ELISA 1:5000-20000 IF 1:50-200
<b>Formulation</b>	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	602
<b>Gene Symbol</b>	BCL3
<b>Uniprot ID</b>	BCL3_HUMAN
<b>Immunogen</b>	Synthesized peptide derived from the human protein at the amino acid range 340-420
<b>Immunogen Region</b>	340-420 aa
<b>Specificity</b>	BCL3 Polyclonal Antibody detects endogenous levels of protein.
<b>Immunogen Sequence</b>	