

Anti-CYP4F8 antibody (80-160 aa) (STJ190288)

STJ190288

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cytochrome P450 4F8 (80-160 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11283
Gene Symbol	CYP4F8
Uniprot ID	CP4F8_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 80-160
Immunogen Region	80-160 aa
Specificity	CP4F8 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	