

Anti-ACOX3 antibody (480-560 aa) (STJ190505)

STJ190505

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Peroxisomal acyl-coenzyme A oxidase 3 (480-560 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8310
Gene Symbol	ACOX3
Uniprot ID	ACOX3_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 480-560
Immunogen Region	480-560 aa
Specificity	ACOX3 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	



Pdcd4 promotes lipid deposition by attenuating PPAR α -mediated fatty acid oxidation in hepatocytes
Mol Cell Endocrinol. 2022 Apr;545:111562. WB Rat 1:1000 liver

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081