

Anti-ADAM30 antibody (400-480 aa) (STJ190513)

STJ190513

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Disintegrin and metalloproteinase domain-containing protein 30 (400-480 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	11085
Gene Symbol	ADAM30
Uniprot ID	ADA30_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 400-480
Immunogen Region	400-480 aa
Specificity	ADA30 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	