

Anti-ADAMTS14 antibody (190-270 aa) (STJ190518)

STJ190518

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-A disintegrin and metalloproteinase with thrombospondin motifs 14 (190-270 aa) is suitable for use in Western Blot and ELISA research applications. |
| Applications | WB/ELISA |
| Host/Source | Rabbit |
| Reactivity | Human/Rat/Mouse |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range | WB 1:500-2000 ELISA 1:5000-20000 |
| Formulation | Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 140766 |
| Gene Symbol | ADAMTS14 |
| Uniprot ID | ATS14_HUMAN |
| Immunogen | Synthesized peptide derived from the human protein at the amino acid range 190-270 |
| Immunogen Region | 190-270 aa |
| Specificity | ATS14 Polyclonal Antibody detects endogenous levels of protein. |
| Immunogen Sequence | |