

Anti-AGXT2 antibody (340-420 aa) (STJ190527)

STJ190527

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Alanine--glyoxylate aminotransferase 2, mitochondrial-3-amino-2-methylpropionate--pyruvate transaminase (340-420 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	64902
Gene Symbol	AGXT2
Uniprot ID	AGT2_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 340-420
Immunogen Region	340-420 aa
Specificity	AGT2 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	