

Anti-ALDH1A3 antibody (290-370 aa) (STJ190531)

STJ190531

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Retinaldehyde dehydrogenase 3 (290-370 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	220
Gene Symbol	ALDH1A3
Uniprot ID	AL1A3_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 290-370
Immunogen Region	290-370 aa
Specificity	AL1A3 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	