

Anti-SCNN1A antibody (320-400 aa) (STJ190541)

STJ190541

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Amiloride-sensitive sodium channel subunit alpha channel subunit alpha (320-400 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6337
Gene Symbol	SCNN1A
Uniprot ID	SCNNA_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 320-400
Immunogen Region	320-400 aa
Specificity	SCNNA Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	