

Anti-ABCB8 antibody (510-590 aa) (STJ190581)

STJ190581

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Mitochondrial potassium channel ATP-binding subunit (510-590 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	11194
Gene Symbol	ABCB8
Uniprot ID	MITOS_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 510-590
Immunogen Region	510-590 aa
Specificity	ABCB8 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	