

Anti-CHD5 antibody (25-75 aa) (STJ190696)

STJ190696

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Chromodomain-helicase-DNA-binding protein 5 (25-75 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	26038
Gene Symbol	CHD5
Uniprot ID	CHD5_HUMAN
Immunogen	Synthesized peptide derived from part of the human protein at the amino acid range 25-75
Immunogen Region	25-75 aa
Specificity	CHD5 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	