

Anti-HECW1 antibody (1-67 aa) (STJ190783)

STJ190783

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-E3 ubiquitin-protein ligase HECW1 (1-67 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23072
Gene Symbol	HECW1
Uniprot ID	HECW1_HUMAN
Immunogen	Synthesized peptide derived from part of the human protein at the amino acid range 1-67
Immunogen Region	1-67 aa
Specificity	HECW1 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	