

Anti-LRP1B antibody (200-280 aa) (STJ190935)

STJ190935

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Low-density lipoprotein receptor-related protein 1B (200-280 aa) is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	IHC-P 1:50-300 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	53353
Gene Symbol	LRP1B
Uniprot ID	LRP1B_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 200-280
Immunogen Region	200-280 aa
Specificity	LRP1B Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	