

Anti-CDHR5 antibody (80-160 aa) (STJ190981)

STJ190981

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cadherin-related family member 5 (80-160 aa) is suitable for use in Western Blot, Immunofluorescence and ELISA research applications.
Applications	WB/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 IF ICC 1:100-500 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	53841
Gene Symbol	CDHR5
Uniprot ID	CDHR5_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 80-160
Immunogen Region	80-160 aa
Specificity	CDHR5 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	