

Anti-NR2C2 antibody (410-490 aa) (STJ191123)

STJ191123

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Nuclear receptor subfamily 2 group C member 2 (410-490 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7182
Gene Symbol	NR2C2
Uniprot ID	NR2C2_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 410-490
Immunogen Region	410-490 aa
Specificity	NR2C2 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	



Western Blot analysis of mouse-brain rat-brain Hela 293 PC-3 cell lysis. primary Antibody was diluted at 1:1000. Secondary antibody (catalog#:STJ99666) was diluted at 1:10000