

Anti-PAIP1 antibody (100-180 aa) (STJ191176)

STJ191176

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Polyadenylate-binding protein-interacting protein 1-binding protein-interacting protein 1 (100-180 aa)
Applications	is suitable for use in Western Blot and ELISA research applications.
Host/Source	WB/ELISA
Reactivity	Rabbit
	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	10605
Gene Symbol	PAIP1
Uniprot ID	PAIP1_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 100-180
Immunogen Region	100-180 aa
Specificity	PAIP1 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	