

Anti-PPP2R5C antibody (360-440 aa) (STJ191383)

STJ191383

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Serine/threonine-protein phosphatase 2A 56 kDa regulatory subunit gamma isoform (360-440 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5527
Gene Symbol	PPP2R5C
Uniprot ID	2A5G_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 360-440
Immunogen Region	360-440 aa
Specificity	2A5G Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	