

Anti-SLC8A2 antibody (380-460 aa) (STJ191438)

STJ191438

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Sodium/calcium exchanger 2/Ca2+-exchange protein 2 (380-460 aa) is suitable for use in Western Blot and ELISA research applications. |
| Applications | WB/ELISA |
| Host/Source | Rabbit |
| Reactivity | Human/Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range | WB 1:500-2000 ELISA 1:5000-20000 |
| Formulation | Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 6543 |
| Gene Symbol | SLC8A2 |
| Uniprot ID | NAC2_HUMAN |
| Immunogen | Synthesized peptide derived from the human protein at the amino acid range 380-460 |
| Immunogen Region | 380-460 aa |
| Specificity | NAC2 Polyclonal Antibody detects endogenous levels of protein. |
| Immunogen Sequence | |