

## Anti-PPAN antibody (STJ191482)

STJ191482

### GENERAL INFORMATION

|                          |                                                                                                                                   |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Suppressor of SWI4 1 homolog is suitable for use in Western Blot and ELISA research applications. |
| <b>Applications</b>      | WB/ELISA                                                                                                                          |
| <b>Host/Source</b>       | Rabbit                                                                                                                            |
| <b>Reactivity</b>        | Human/Mouse                                                                                                                       |

### PRODUCT PROPERTIES

|                            |                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                                            |
| <b>Clone ID</b>            |                                                                                                                       |
| <b>Concentration</b>       | 1 mg/mL                                                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                                                          |
| <b>Purification</b>        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Dilution Range</b>      | WB 1:500-2000<br>ELISA 1:5000-20000                                                                                   |
| <b>Formulation</b>         | Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.                                                         |
| <b>Isotype</b>             | IgG                                                                                                                   |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

### TARGET INFORMATION

|                           |                                                                |
|---------------------------|----------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">56342/692312</a>                                   |
| <b>Gene Symbol</b>        | <a href="#">PPAN</a>                                           |
| <b>Uniprot ID</b>         | <a href="#">SSF1_HUMAN</a>                                     |
| <b>Immunogen</b>          | Synthesized peptide derived from part of the human protein     |
| <b>Immunogen Region</b>   |                                                                |
| <b>Specificity</b>        | SSF1 Polyclonal Antibody detects endogenous levels of protein. |
| <b>Immunogen Sequence</b> |                                                                |

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081