

Anti-SPAG9 antibody (590-670 aa) (STJ191696)

STJ191696

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-C-Jun-amino-terminal kinase-interacting protein 4 (590-670 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9043
Gene Symbol	SPAG9
Uniprot ID	JIP4_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 590-670
Immunogen Region	590-670 aa
Specificity	JIP4 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	