

Anti-GRIP2 antibody (780-860 aa) (STJ191851)

STJ191851

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Glutamate receptor-interacting protein 2 (780-860 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	80852
Gene Symbol	GRIP2
Uniprot ID	GRIP2_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 780-860
Immunogen Region	780-860 aa
Specificity	GRIP2 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	