

Anti-SMARCD1 antibody (268-318 aa) (STJ191929)

STJ191929

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily D member 1 (268-318 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Range	ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6602
Gene Symbol	SMARCD1
Uniprot ID	SMRD1_HUMAN
Immunogen	Synthesized peptide derived from part of the human protein at the amino acid range 268-318
Immunogen Region	268-318 aa
Specificity	SMRD1 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	