

Anti-CYP4A22 antibody (403-453 aa) (STJ191942)

STJ191942

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cytochrome P450 4A22 (403-453 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	284541
Gene Symbol	CYP4A22
Uniprot ID	CP4AM_HUMAN
Immunogen	Synthesized peptide derived from part of the human protein at the amino acid range 403-453
Immunogen Region	403-453 aa
Specificity	CP4AM Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	