

## Anti-LRIG3 antibody (STJ192075)

STJ192075

### GENERAL INFORMATION

|                          |                                                                                                                                                                     |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                                                  |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Leucine-rich repeats and immunoglobulin-like domains protein 3 is suitable for use in Western Blot and ELISA research applications. |
| <b>Applications</b>      | WB/ELISA                                                                                                                                                            |
| <b>Host/Source</b>       | Rabbit                                                                                                                                                              |
| <b>Reactivity</b>        | Human/Mouse                                                                                                                                                         |

### PRODUCT PROPERTIES

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>      | Polyclonal                                                                                                            |
| <b>Clone ID</b>       |                                                                                                                       |
| <b>Concentration</b>  | 1 mg/mL                                                                                                               |
| <b>Conjugation</b>    | Unconjugated                                                                                                          |
| <b>Purification</b>   | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Dilution Range</b> | WB 1:500-2000<br>ELISA 1:5000-20000                                                                                   |
| <b>Formulation</b>    | Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.                                                         |
| <b>Isotype</b>        | IgG                                                                                                                   |
| <b>Storage</b>        | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |
| <b>Instruction</b>    |                                                                                                                       |

### TARGET INFORMATION

|                           |                                                                 |
|---------------------------|-----------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">121227</a>                                          |
| <b>Gene Symbol</b>        | <a href="#">LRIG3</a>                                           |
| <b>Uniprot ID</b>         | <a href="#">LRIG3_HUMAN</a>                                     |
| <b>Immunogen</b>          | Synthesized peptide derived from part of the human protein      |
| <b>Immunogen Region</b>   |                                                                 |
| <b>Specificity</b>        | LRIG3 Polyclonal Antibody detects endogenous levels of protein. |
| <b>Immunogen Sequence</b> |                                                                 |