

Anti-PLPPR4 antibody (55-80 aa) (STJ192133)

STJ192133

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospholipid phosphatase-related protein type 4 (55-80 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9890
Gene Symbol	PLPPR4
Uniprot ID	PLPR4_HUMAN
Immunogen	Synthesized peptide derived from part of the human protein.the amino acid55-80
Immunogen Region	55-80 aa
Specificity	LPFR4 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	