

Anti-TNFRSF21 antibody (430-510 aa) (STJ192397)

STJ192397

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Tumor necrosis factor receptor superfamily member 21 (430-510 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	27242
Gene Symbol	TNFRSF21
Uniprot ID	TNR21_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 430-510
Immunogen Region	430-510 aa
Specificity	TNR21 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	