

Anti-IGSF3 antibody (320-400 aa) (STJ192411)

STJ192411

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Immunoglobulin superfamily member 3 (320-400 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3321
Gene Symbol	IGSF3
Uniprot ID	IGSF3_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 320-400
Immunogen Region	320-400 aa
Specificity	IGSF3 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	