

Anti-ZC3HAV1 antibody (211-260 aa) (STJ192605)

STJ192605

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Zinc finger CCCH-type antiviral protein 1 (211-260 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	56829
Gene Symbol	ZC3HAV1
Uniprot ID	ZCCHV_HUMAN
Immunogen	Synthesized peptide derived from the human protein at the amino acid range 211-260
Immunogen Region	211-260 aa
Specificity	ZCCHV Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	