

Anti-CLTC antibody (482-532 aa) (STJ192903)

STJ192903

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Clathrin heavy chain 1 (482-532 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1213
Gene Symbol	CLTC
Uniprot ID	CLH1_HUMAN
Immunogen	Synthesized peptide derived from part of the human protein at the amino acid range 482-532
Immunogen Region	482-532 aa
Specificity	CLH1 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	