

Anti-ANKS1A antibody (44-94 aa) (STJ192929)

STJ192929

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ankyrin repeat and SAM domain-containing protein 1A (44-94 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	23294
Gene Symbol	ANKS1A
Uniprot ID	ANS1A_HUMAN
Immunogen	Synthesized peptide derived from part of the human protein at the amino acid range 44-94
Immunogen Region	44-94 aa
Specificity	ANS1A Polyclonal Antibody detects endogenous levels of protein.
Immunogen Sequence	