

Anti-Phospho-PDGFRB-Tyr857 antibody (STJ193215)

STJ193215

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Platelet-derived growth factor receptor beta-Tyr857 is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	IHC/IF/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-P 1:50-200
Range	WB 1:500-2000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5159
Gene Symbol	PDGFRB
Uniprot ID	PDGFRB_HUMAN
Immunogen	Synthesized peptide derived from the human PDGFRb (Phospho-Tyr857)
Immunogen Region	
Specificity	This antibody detects endogenous phospho levels of PDGFRb (Phospho-Tyr857) at Human:Y857, Mouse:Y856, Rat:Y856
Immunogen Sequence	