

Anti-SARS-nCoV-2 antibody [S03FC] (STJ193235)

STJ193235

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-SARS-nCov-2 Spike is suitable for use in ELISA, ECL and Immunogold research applications.
Applications	ELISA/ECL/Immunogold
Host/Source	Human
Reactivity	COVID-19

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S03FC
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from 293T expression by affinity-chromatography using specific immunogen.
Dilution Range	ELISA 1:5000-1:100000 WB 1:1000-1:10000
Formulation	Liquid in PBS
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	The recombinant SARS-CoV-2 (2019-nCoV) Spike Protein S
Immunogen Region	
Specificity	This antibody detects SARS-CoV-2 (2019-nCoV) Spike Protein S1 domain
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081