

## Anti-CIPC antibody (141-191 aa) (STJ194483)

STJ194483

### GENERAL INFORMATION


|                          |                                                                                                                                     |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                  |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-CLOCK-interacting pacemaker (141-191 aa) is suitable for use in Western Blot research applications. |
| <b>Applications</b>      | WB                                                                                                                                  |
| <b>Host/Source</b>       | Rabbit                                                                                                                              |
| <b>Reactivity</b>        | Human/Mouse                                                                                                                         |

### PRODUCT PROPERTIES

|                            |                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                                            |
| <b>Clone ID</b>            |                                                                                                                       |
| <b>Concentration</b>       | 1 mg/mL                                                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                                                          |
| <b>Purification</b>        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Dilution Range</b>      | WB 1:500-2000                                                                                                         |
| <b>Formulation</b>         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                               |
| <b>Isotype</b>             | IgG                                                                                                                   |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

### TARGET INFORMATION

|                           |                                                                                  |
|---------------------------|----------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 85457                                                                            |
| <b>Gene Symbol</b>        | CIPC                                                                             |
| <b>Uniprot ID</b>         | CIPC_HUMAN                                                                       |
| <b>Immunogen</b>          | Synthesized peptide derived from the human K1737 at the amino acid range 141-191 |
| <b>Immunogen Region</b>   | 141-191 aa                                                                       |
| <b>Specificity</b>        | This antibody detects endogenous levels of K1737 at Human/Mouse                  |
| <b>Immunogen Sequence</b> |                                                                                  |



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4A°C over night

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081